

Instructions for Use – Screening Cells in CellStab



IVD

REAGENT RED CELLS *for the detection of red cell antibodies*

For in vitro diagnostic use only
Product codes PR102 and PR122
NBS Reagents, NHSBT Liverpool,
14 Estuary Banks, Liverpool L24 8RB
Tel 0151 268 7157



Blood and Transplant

Intended use

These reagent red blood cells are intended to be used to screen plasma/serum samples for presence of red cell antibodies by serological means.

Principles of the examination method

Plasma/serum samples are incubated with reagent red cells to determine the presence of agglutinins by direct and/or indirect methods

Components

These cells are supplied as a 0.8±0.2% suspension to be used directly from the vial with BioRad ID cards, in accordance with the manufacturer's package insert. These reagent red cells, prepared from non-remunerated blood donors, are leucodepleted, washed and suspended in a preservative solution – CellStab solution, which has been specially formulated to preserve red cell integrity and antigenicity, consisting of glycine buffered saline, containing sugars, trimethoprim and sulfamethoxazol as preservatives.

Reagent Preparation

Mix before use

Storage and shelf life after first opening

Store at 2°- 8°C

Do not freeze

Do not use beyond the notified expiry date

Warnings and precautions

For professional use only.

For use only with BioRad -ID gel cards.

It is imperative to use accurate, properly calibrated volumetric pipettes in BioRad ID systems to avoid variations which may affect the test outcome.

For red cells that have been treated with the enzyme papain the following antigens will be absent or reduced: M, N, S, s Fya, Fyb, Ch/Rg, In, JMH, Xga.

Some loss of antigenic expression may occur during the stated shelf life. Since this loss cannot be predicted or controlled and is partly determined by the characteristics of individual blood donations or donors, the recommended conditions of storage and use must be rigidly applied.

Do not use if red cells appear obviously discoloured or haemolysed.

Cells must not be pooled

The donations used in this product have been tested at source and found negative for the mandatory microbiological tests required by the UK BTS at the time of donation. No known test methods can offer assurances that products derived from human blood will not transmit infectious diseases. Appropriate care should be taken in the use and disposal of this product.

Primary sample collection, handling and storage

Clotted serum or EDTA plasma samples may be used, usually less than 7 days old, different time scales apply to recently transfused patients

Examination procedure

It is important to detect the presence of clinically significant red cell antibodies in a patients serum/plasma and subsequently identify the specificity of such antibodies in order to ensure that any subsequent transfusion is as free from risk of a red cell transfusion reaction as possible within the limits of the techniques used. Incubating the patients serum/plasma and NBS Reagents antibody screening cells in accordance with the BioRad package insert assists in this aim.

Refer to the package insert of the cards being used; only this technique should be used. NBS Reagents antibody screening cells can be used in place of the ID-DiaCell reagents recommended in the package insert.

Control procedure

Each batch of tests should be controlled with suitable positive and negative controls

Interpretation of results

The presence of agglutination indicates a positive result, meaning that an antibody may be present in the sample and may only be seen against 1 cell of the set and will require further investigation to assign specificity to the reaction seen.

Performance characteristics

The reagent red cells selected to be used for the detection of antibodies are negative for Wra, and positive for Lub and Kpb unless stated.

The antigenic status of these red cells has been determined using, wherever possible, at least two examples of antisera directed against that antigen. The designation of positive or negative status for a particular antigen relates to the normal expression of that antigen, if an individual cell is known to possess a weak or variant form of an antigen, this is indicated on the profile

Limitations of the examination procedure

If controls set up with the batch of tests fail to give required results then all tests must be repeated.

Literature references

These reagents comply with

- The requirements of Directive 98/79/EC on in vitro diagnostic medical devices
- The recommendations contained in current version of the Guidelines for the Blood Transfusion Services in the UK.
- BCSH Guidelines for compatibility procedures - current version.