

## Instructions for Use - rr screening cells



Reagent red cells for the  
detection of red cell antibodies  
in patients who have received  
prophylactic Anti-D

For *in vitro* diagnostic use only

Product codes: PR106 and  
PR107

NHSBT Reagents,  
14 Estuary Banks, Liverpool  
L24 8RB

Tel 0151 268 7157



Blood and Transplant

### Intended use

These reagent red blood cells are intended to be used to screen plasma/serum samples for presence of red cell antibodies by serological means.

### Principles of the examination method

Plasma/serum samples are incubated with reagent red cells to determine the presence of agglutinins by direct and/or indirect methods

### Components

These cells are supplied as a  $0.8 \pm 0.2\%$  suspension to be used directly from the vial with BioRad cards, in accordance with the manufacturer's package insert.

These reagent red cells, prepared from non-remunerated blood donors, are leucodepleted, washed and suspended in a preservative solution – CellStab which has been specially formulated to preserve red cell integrity and antigenicity, consisting of glycine buffered saline, containing sugars, trimethoprim and sulfamethoxazol as preservatives.

### Reagent Preparation

Mix before use.

### Storage and shelf life after first opening

Store at  $2^{\circ}$ -  $8^{\circ}$ C.

Do not freeze.

Do not use beyond the notified expiry date.

### Warnings and precautions

For professional use only.

This screening set must not be used for general antibody screening.

This screening cell set is intended for use with patients who have received prophylactic Anti-D as part of their antenatal care

For use only with BioRad gel cards.

It is imperative to use accurate, properly calibrated volumetric pipettes in BioRad systems to avoid variations which may affect the test outcome.

Some loss of antigenic expression may occur during the stated shelf life. Since this loss cannot be predicted or controlled and is partly determined by the characteristics of individual blood donations or donors, the recommended conditions of storage and use must be rigidly applied.

Do not use if red cells appear obviously discoloured or haemolysed.

Cells must not be pooled.

The donations used in this product have been tested at source and found negative for the mandatory microbiological tests required by the Guidelines for UK BTS at the time of donation. No known test methods can offer assurances that products derived from human blood will not transmit infectious diseases. Appropriate care should be taken in the use and disposal of this product.

### Primary sample collection, handling and storage

Clotted serum or EDTA plasma samples may be used, usually less than 7 days old, different time scales apply to recently transfused patients

### Examination procedure

It is important to detect the presence of clinically significant red cell antibodies in a patients serum/plasma and subsequently identify the specificity of such antibodies in order to ensure that any subsequent transfusion is as free from risk of a red cell transfusion reaction as possible within the limits of the techniques used. Incubating the patients serum/plasma and Reagents antibody screening cells in accordance with the BioRad ID package insert assists in this aim.

### Control procedure

Each batch of tests should be controlled with suitable positive and negative controls

### Interpretation of results

The presence of agglutination indicates a positive result, meaning that an antibody may be present in the sample and may only be seen against 1 cell of the set and will require further investigation to assign specificity to the reaction seen.

### Performance characteristics

The rr (PR106) reagent red cells selected to be used for the detection of antibodies are negative for Wra, and positive for Lub and Kpb unless stated.

The r'r and r'r subset (PR107) reagent red cells selected to be used for the detection of antibodies are negative for Wra with one of them being positive for Kpa unless otherwise stated.

The antigenic status of these red cells has been determined using, wherever possible, at least two examples of antisera directed against that antigen.

The designation of positive or negative status for a particular antigen relates to the normal expression of that antigen, if an individual cell is known to possess a weak or variant form of an antigen, this is indicated on the profile

### Limitations of the examination procedure

If controls set up with the batch of tests fail to give required results then all tests must be repeated.

### Literature references

These reagents comply with:

The requirements of Directive 98/79/EC on in vitro diagnostic medical devices

The recommendations contained in current version of the Guidelines for the Blood Transfusion Services in the UK. BSH Guidelines for compatibility procedures - current version.



Quality First International OÜ, Laki 30, 12915 Tallinn, Estonia