

Note

Transportation of Blood Components

Container capacities include up to six neo natal units as one adult unit for all component types.

Fill any excess air space between PCMs/Dry ice and lid and around units with paper towel/bubble wrap or paper towel/bubble wrap cushions to ensure no movement of components during transportation.

All components when packed in Va-Q-tecs must be packed flat to decrease the likelihood of a pack splitting. If packed vertically force would be projected along the seam of a pack making it more likely to split.

Transportation of Blood Donations

Va-Q-Tec transport containers **MUST** only be used for transporting blood donations from stock holding units to manufacturing centres when a bulk movement vehicle unexpectedly becomes unavailable.

Whole blood donations must be laid flat with the satellite bags and filters underneath. Units can overlap in order to fit into the containers.

Fill any excess air space between PCMs and lid and around units with paper towel/bubble wrap or paper towel/bubble wrap cushions to ensure no movement of donations during transportation.

Extreme Weather Protocol

At those times when external temperatures **FALL BELOW -5°C** or **RISE ABOVE +35°C** only, the **EXTREME WEATHER PROTOCOL** listed on this datasheet **MUST** be instigated for transportation of all blood components.

The protocol covers external temperatures down to -10°C or up to +40°C.

DAT48/14 – Capacity and Transportation Time Limits for Transport Containers



Blood and Transplant

Copy No:

Effective date: 28MAY2025

VA-Q-TEC LARGE TRANSPORT CONTAINERS

Component Type	Maximum Capacity (Units)	Maximum Time Units out of Controlled Storage in Container (Hours)	Temperature Stabilisation Material Type	Configuration of Temperature Stabilisation Material	External Temperature
Red Blood Cells	12	9	Blue (+4°C +/-2°C)	1 Bottom 2 Opposite Sides 1 Top	-5°C to +35°C
Platelets	15	8	Green (+22°C +/-2°C)		
Clinical Buffy Coats	10				
Granulocytes	2				

VA-Q-TEC SMALL TRANSPORT CONTAINERS

Component Type	Maximum Capacity (Units)	Maximum Time Units out of Controlled Storage in Container (Hours)	Temperature Stabilisation Material Type	Configuration of Temperature Stabilisation Material	External Temperature
Red Blood Cells	6	5.5	Blue (+4°C +/-2°C)	2 bottom 2 top Side by side	-5°C to +35°C
Platelets		7	Green (+22°C +/-2°C)		
Adult FFP	4	11	Dry Ice	2 x 500g bags Top 1 x 500g bag Bottom	
Low Volume Frozen Components	6				

VA-Q-TEC MEDIUM TRANSPORT CONTAINERS

Component Type	Maximum Capacity (Units)	Maximum Time Units out of Controlled Storage in Container (Hours)	Temperature Stabilisation Material Type	Configuration of Temperature Stabilisation Material	External Temperature
Red Blood Cells	15	3	Blue (+4°C +/-2°C)	2 bottom 2 top Side by side	0°C to +30°C
Platelets		5	Green (+22°C +/-2°C)		
FFP and Low Volume Frozen Components	10	9.5	Dry ice	2 x 500g bags Top 1 x 500g bag Bottom	-5°C to +35°C

Controlled if copy number stated on document and issued by QA

(Template Version 03/02/2020)

Cross-Referenced in Primary Document: SOP4147

DAT48/14 – Capacity and Transportation Time Limits for Transport Containers



Blood and Transplant

Copy No:

Effective date: 28MAY2025

EXTREME WEATHER PROTOCOL (-10°C to +40°C)

This protocol should only be instigated if external temperatures are below -5°C and above +35°C

VA-Q-TEC LARGE TRANSPORT CONTAINERS

Component Type	Maximum Capacity (Units)	Maximum Time Units out of Controlled Storage in Container (Hours)	Temperature Stabilisation Material Type	Configuration of Temperature Stabilisation Material	External Temperature
Red Blood Cells	10*	9	Blue (+4°C +/-2°C)	1 Bottom 2 Opposite Sides 1 Top	-10°C to +40°C
Platelets	15	8	Green (+22°C +/-2°C)		
Note *For the Extreme Weather Protocol only, maximum capacity is 10 units of Red Blood Cells					

VA-Q-TEC SMALL TRANSPORT CONTAINERS

Component Type	Maximum Capacity (Units)	Maximum Time Units out of Controlled Storage in Container (Hours)	Temperature Stabilisation Material Type	Configuration of Temperature Stabilisation Material	External Temperature
Red Blood Cells	6	5.5	Blue (+4°C +/-2°C)	2 bottom 2 top Side by side	-10°C to +40°C
Platelets	6	3	Green (+22°C +/-2°C)		
Frozen Components	4	11	Dry Ice	2 x 500g bags Top **2 x 500g bag Bottom	-10°C to +40°C
<p>Note</p> <p>**For the Extreme Weather Protocol only, two 500g of dry ice must be placed at the bottom of the container</p>					

VA-Q-TEC MEDIUM TRANSPORT CONTAINERS

Component Type	Maximum Capacity (Units)	Maximum Time Units out of Controlled Storage in Container (Hours)	Temperature Stabilisation Material Type	Configuration of Temperature Stabilisation Material	External Temperature
Frozen Components	10	10	Dry Ice	2 x 500g bags Top **2 x 500g bag Bottom	-10°C to +40°C
<p>Note</p> <p>**For the Extreme Weather Protocol only, two 500g of dry ice must be placed at the bottom of the container</p>					







Controlled if copy number stated on document and issued by QA

(Template Version 03/02/2020)








Cross-Referenced in Primary Document: SOP4147

PACKING CONFIGURATIONS

Diagrams are for illustration purposes only and are not to scale


Packing Configuration Examples for SMALL transport containers (Minimum and Maximum load) <i>The configurations shown are only examples. Other configurations are valid if the units are lying flat in stacks equally distributed in the container. The units can overlap if there is limited space.</i>	
1 unit Red Blood Cells – Top View of Container  <p>Ports Folded Laid flat 1 stack x 1 unit</p>	6 units Red Blood Cells – Top View of Container  <p>Ports Folded Laid flat 2 stacks x 3 units</p>
1 unit Platelets - Top View of Container  <p>Laid flat 1 stack x 1 unit</p>	6 units Platelets – Top View of Container  <p>Laid flat 1 stack x 6 units</p>
1 unit FFP – Top View of Container  <p>Laid flat 1 stack x 1 unit</p>	4 units FFP – Top View of Container  <p>Laid flat 1 stack x 4 units</p>

DAT48/14 – Capacity and Transportation Time Limits for Transport Containers

Packing Configuration Examples for MEDIUM transport containers (Minimum and Maximum load) The configurations shown are only examples. Other configurations are valid if the units are lying flat in stacks equally distributed in the container. The units can overlap if there is limited space.	
1 and 7 units Red Blood Cells – Top View of Container  Ports Folded Laid flat 1 stack x 1 unit, and  3 stacks x 2 units 1 stack x 1 unit	15 units Red Blood Cells – Top View of Container  Ports Folded Laid flat 3 stacks x 4 units 1 stack x 3 units
1 unit Platelets – Top View of Container  Laid Flat 1 stack x 1 units	15 units Platelets – Top View of Container  Laid Flat 1 stack x 15 units
1 unit FFP Top View of Container  Laid Flat 1 stack x 1 unit	10 units FFP – Top View of Container  Laid Flat 1 stack x 10 units


Packing Configuration Examples for LARGE TRANSPORT CONTAINERS (Minimum and Maximum load)
The configurations shown are only examples. Other configurations are valid if the units are lying flat in stacks equally distributed in the container. The units can overlap if there is limited space

1 unit Red Blood Cells – Top View of Container




Ports Folded
Laid Flat
1 stack x 1 unit

12 units Red Blood Cells – Top View of Container




Ports Folded
Laid flat
2 stacks x 6 units

1 unit Platelets – Top View of Container



Laid Flat
1 stack x 1 unit

15 units Platelets – Top View of Container



Laid Flat
1 stack x 15 units

Continued Packing Configuration Examples for **SMALL** transport containers

The configurations shown are only examples. Other configurations are valid if the units are lying flat in stacks equally distributed in the container. The units can overlap if there is limited space

6 units of red blood cells

Ports Folded
Laid flat
2 stack x 3 unit
With overlap



6 units of Platelets

Laid Flat
1 stack x 6 units



4 units of Plasma
Laid flat
1 stack x 4 units



Continued Packing Configuration Examples for MEDIUM transport containers

The configurations shown are only examples. Other configurations are valid if the units are lying flat in stacks equally distributed in the container. The units can overlap if there is limited space

15 units of red blood cells

Ports Folded
Laid flat
3 stack x 4 units
1 Stack x 3 units



8 units of Platelets

Laid Flat
1 stack x 8 units



10 units of Plasma
Laid flat
2 stack x 5 units
With overlap



DAT48/14 – Capacity and Transportation Time Limits for Transport Containers

Continued Packing Configuration Examples for LARGE transport containers

The configurations shown are only examples. Other configurations are valid if the units are lying flat in stacks equally distributed in the container. The units can overlap if there is limited space.

12 units of red
blood cells

Ports Folded
Laid flat
4 stacks x 3 units
With overlap



8 units of Platelets

Laid Flat
1 stack x 8 units



TRANSPORTATION OF BLOOD DONATIONS

VA-Q-TEC LARGE TRANSPORT CONTAINERS

Donation Type	Maximum Capacity (Donations)	Maximum Time Units out of Controlled Storage in Container (Hours)	Temperature Stabilisation Material Type	Configuration of Temperature Stabilisation Material	External Temperature
CD Platelets One Donation (2 split packs)	7	8	Green (+22°C +/-2°C)	1 Bottom 2 Opposite Sides 1 Top	-5°C to +35°C
Donated Plasma					
Whole Blood Donation	4				

VA-Q-TEC SMALL TRANSPORT CONTAINERS

Donation Type	Maximum Capacity (Donations)	Maximum Time Units out of Controlled Storage in Container (Hours)	Temperature Stabilisation Material Type	Configuration of Temperature Stabilisation Material	External Temperature
CD Platelets One Donation (2 split packs)	3	7	Green (+22°C +/-2°C)	2 Bottom 2 Top Side by Side	-5°C to +35°C
Donated Plasma					
Whole Blood Donation		5			

VA-Q-TEC MEDIUM TRANSPORT CONTAINERS

Donation Type	Maximum Capacity (Donations)	Maximum Time Units out of Controlled Storage in Container (Hours)	Temperature Stabilisation Material Type	Configuration of Temperature Stabilisation Material	External Temperature
CD Platelets One Donation (2 split packs)	7	5	Green (+22°C +/-2°C)	2 Bottom 2 Top Side by Side	0°C to +30°C
Donated Plasma					
Whole Blood Donation	9	3.5			