

ID Panel Profile



0843
IVD



Blood and Transplant

NHSBT REAGENTS Panel 1

Product	Lot No.	Product	Lot No.	Product	Lot No.
ID Panel in Alsevers	R144 3434	ID Panel Papainised in Alsevers	R154 3434	ID Panel in LISP	R146 3434
ID Panel in CellStab	R143 3434	ID Panel Papainised in CellStab	R153 3434	Expiry Date : 2019.08.15	
ID Panel in CellMedia	R163 3434	ID Panel Papainised in CellMedia	R173 3434		

Patient's Name		Ref. No.	Sample No.	Conclusion	
D.O.B.				Tested by	Date

Unless otherwise indicated, all cells are positive for Kp^b and Lu^b and negative for Wr^a and Co^b.

Instructions for use can be found at <http://www.blood.co.uk/reagents>

	Rh	C	D	E	c	e	C ^w	M	N	S	s	P1	Lu ^a	K	k	Kp ^a	Le ^a	Le ^b	Fy ^a	Fy ^b	Jk ^a	Jk ^b	Other																
1	R ₁ ^w R ₁	+	+	0	0	+	+	0	+	0	+	0	+	0	+	0	0	+	+	0	+	0																	
2	R ₁ R ₁	+	+	0	0	+	0	+	0	+	0	1	0	+	+	0	0	+	0	+	0	+																	
3	R ₂ R ₂	0	+	+	+	0	0	0	+	0	+	0	0	0	+	0	0	0	+	0	0	+																	
4	r'r	+	0	0	+	+	0	+	0	0	+	2	0	0	+	0	+	0	+	0	0	+																	
5	r''r	0	0	+	+	+	0	+	0	+	0	0	0	0	+	0	0	+	0	+	+	0																	
6	rr	0	0	0	+	+	0	+	+	0	+	4	0	+	0	0	0	+	+	0	+	0																	
7	rr	0	0	0	+	+	0	0	+	0	+	1	0	+	+	0	+	0	0	+	+	0																	
8	rr	0	0	0	+	+	0	+	0	+	0	0	0	0	+	+	0	0	0	+	0	+																	
9	rr	0	0	0	+	+	0	0	+	0	+	0	+	0	+	0	+	0	0	+	+	0																	
10	rr	0	0	0	+	+	0	0	+	+	0	3	0	0	+	0	0	0	+	0	0	+																	
Reagent																																							
Lot No.																																							