

3 Cell Screen Product Profile



Blood and Transplant

IVD

NHSBT REAGENTS

Product	Lot No	Product	Lot No.	Product	Lot No.	Expiry Date
Alsevers	R121 3565	CellStab	R122 3565	CellMedia	R123 3565	2019.08.22

Unless otherwise indicated, all cells are positive for Kp^b and Lu^b and negative for Wr^a, Lu^a and Co^b

Instructions for use can be found at www.blood.co.uk/reagents

	Rh	C	D	E	c	e	C ^w	M	N	S	s	P1	K	k	Kp ^a	Le ^a	Le ^b	Fy ^a	Fy ^b	Jk ^a	Jk ^b	Other
1	R₁^wR₁	+	+	0	0	+	+	+	0	+	0	+	0	+	0	0	+	+	+	0	+	
2	R₂R₂	0	+	+	+	0	0	0	+	0	+	+	+	+	0	+	0	+	0	0	+	
3	rr	0	0	0	+	+	0	+	+	0	+	+	0	+	+	0	+	0	+	+	0	