

2 Cell Screen Profile PR101, PR102 & PR103



NHSBT REAGENTS

Product	Lot No.	Product	Lot No.
Alsevers	R101 3479	CellStab	R102 3479
CellMedia	R103 3479	EXPIRY DATE: 2019.08.29	

Unless otherwise indicated, all cells are positive for Kp^b and Lu^b and negative for Wr^a , Lu^a and Co^b
 Instructions for use can be found at www.blood.co.uk/reagents

	Rh	C	D	E	c	e	C ^w	M	N	S	s	P1	K	k	Kp ^a	Le ^a	Le ^b	Fy ^a	Fy ^b	Jk ^a	Jk ^b	Other
1	R ₁ R ₁	+	+	0	0	+	0	+	+	+	0	+	+	+	0	0	+	0	+	+	0	Yka- Ytb+
2	R ₂ R ₂	0	+	+	+	0	0	+	0	0	+	+	0	+	0	+	0	+	0	0	+	Cob+