

2 Cell Screen Profile PR101, PR102 & PR103



NHSBT REAGENTS

Product	Lot No.	Product	Lot No.
Alsevers	R101 3480	CellStab	R102 3480
CellMedia	R103 3480	EXPIRY DATE: 2019.09.26	

Unless otherwise indicated, all cells are positive for Kp^b and Lu^b and negative for Wr^a , Lu^a and Co^b
 Instructions for use can be found at www.blood.co.uk/reagents

	Rh	C	D	E	c	e	C ^w	M	N	S	s	P1	K	k	Kp ^a	Le ^a	Le ^b	Fy ^a	Fy ^b	Jk ^a	Jk ^b	Other
1	R₁R₁	+	+	0	0	+	0	+	0	0	+	+	+	+	0	+	0	0	+	+	0	
2	R₂R₂	0	+	+	+	0	0	+	+	+	0	+	0	+	0	0	+	+	0	0	+	