

ID Panel Profile



Blood and Transplant

NHSBT REAGENTS Panel 1

Product	Lot No.	Product	Lot No.	Product	Lot No.
ID Panel in Alsevers	R144 3462	ID Panel Papainised in Alsevers	R154 3462	ID Panel in LISP	R146 3462
ID Panel in CellStab	R143 3462	ID Panel Papainised in CellStab	R153 3462	Expiry Date : 2020.09.24	
ID Panel in CellMedia	R163 3462	ID Panel Papainised in CellMedia	R173 3462		

Patient's Name		Ref. No.	Sample No.	Conclusion
D.O.B.				Tested by Date

Unless otherwise indicated, all cells are positive for Kp^b and Lu^b and negative for Wr^a and Co^b.

Instructions for use can be found at <http://www.blood.co.uk/reagents>

EC REP Quality First International OÜ, Laki 30, 12915 Tallinn, Estonia

	Rh	C	D	E	c	e	C ^w	M	N	S	s	P1	Lu ^a	K	k	Kp ^a	Le ^a	Le ^b	Fy ^a	Fy ^b	Jk ^a	Jk ^b	Other						
1	R ₁ ^w R ₁	+	+	0	0	+	+	0	+	0	+	0	0	0	+	0	0	+	0	+	0	+							
2	R ₁ R ₁	+	+	0	0	+	0	+	0	+	0	3	0	+	+	0	+	0	+	0	+	0							
3	R ₂ R ₂	0	+	+	+	0	0	0	+	0	+	4	0	0	+	0	0	+	+	0	0	+							
4	r'r	+	0	0	+	+	0	+	0	0	+	2	0	0	+	0	+	0	+	0	0	+							
5	r''r	0	0	+	+	+	0	+	0	+	0	3	0	0	+	0	+	0	0	+	+	0							
6	rr	0	0	0	+	+	0	+	+	0	+	0	0	+	0	0	0	+	0	w	0	+							
7	rr	0	0	0	+	+	0	0	+	0	+	1	0	+	+	0	+	0	0	+	+	0							
8	rr	0	0	0	+	+	0	0	+	0	+	0	0	0	+	+	0	+	+	0	+	0							
9	rr	0	0	0	+	+	0	0	+	+	0	3	0	0	+	0	0	+	+	0	0	+							
10	rr	0	0	0	+	+	0	+	0	0	+	4	+	0	+	0	0	+	+	0	+	0	Rg-						
Reagent																	DAT Profile	Anti-IgG	IgA	IgM	C3c	C3d	Ctrl						
Lot No.																	Result												