

rr Screen Profile Products PR106, PR107, PR108 & PR109



Blood and Transplant

NHSBT REAGENTS

Product	Lot No.	Expiry
CellStab	R106 3230 R107 3230	2020.11.05
CellMedia	R108 3230 R109 3230	

Unless otherwise indicated, all cells are positive for Kp^b and Lu^b and negative for Wr^a, Lu^a and Co^b

Instructions for use can be found at www.blood.co.uk/reagents

		Rh	C	D	E	c	e	C ^w	M	N	S	s	P1	K	k	Kp ^a	Le ^a	Le ^b	Fy ^a	Fy ^b	Jk ^a	Jk ^b	Other
PR106/ PR108	1	rr	0	0	0	+	+	0	+	0	+	0	0	+	+	0	+	0	+	0	0	+	
	2	rr	0	0	0	+	+	0	0	+	0	+	+	0	+	0	0	+	0	+	+	0	Cob+
PR107/ PR109	3	r'r	+	0	0	+	+	0	0	+	0	+	+	0	+	+	+	0	0	+	+	0	
	4	r''r	0	0	+	+	+	0	+	+	+	0	0	0	+	0	0	+	+	0	0	+	

EC REP Quality First International OÜ, Laki 30, 12915 Tallinn, Estonia

(Template Version 01/04)