

**rr Screen Profile Products PR106, PR107, PR108 & PR109****Blood and Transplant****NHSBT REAGENTS**

Product	Lot No.	Expiry
<b>CellStab</b>	<b>R106 3241</b> <b>R107 3241</b>	<b>2021.09.09</b>
<b>CellMedia</b>	<b>R108 3241</b> <b>R109 3241</b>	

**Unless otherwise indicated, all cells are positive for Kp<sup>b</sup> and Lu<sup>b</sup> and negative for Wr<sup>a</sup>, Lu<sup>a</sup> and Co<sup>b</sup>**

Instructions for use can be found at [www.blood.co.uk/reagents](http://www.blood.co.uk/reagents)

		Rh	C	D	E	c	e	C <sup>w</sup>	M	N	S	s	P1	K	k	Kp <sup>a</sup>	Le <sup>a</sup>	Le <sup>b</sup>	Fy <sup>a</sup>	Fy <sup>b</sup>	Jk <sup>a</sup>	Jk <sup>b</sup>	Other
PR106/ PR108	1	rr	0	0	0	+	+	0	+	0	+	0	0	+	+	0	+	0	+	0	0	+	
	2	rr	0	0	0	+	+	0	0	+	0	+	+	0	+	0	0	+	0	+	+	0	
PR107/ PR109	3	r'r	+	0	0	+	+	0	0	+	0	+	+	0	+	0	0	+	0	+	+	0	Cob+
	4	r''r	0	0	+	+	+	0	+	0	+	0	+	+	+	+	+	0	+	w	+	+	

**EC REP** Quality First International OÜ, Laki 30, 12915 Tallinn, Estonia

(Template Version 01/04)