

3 Cell Screen Product Profile



NHSBT REAGENTS

Product	Lot No	Product	Lot No.	Product	Lot No.	Expiry Date
Alsevers	R121 3693	CellStab	R122 3693	CellMedia	R123 3693	2024.06.20

Unless otherwise indicated, all cells are positive for Kp^b and Lu^b and negative for Wr^a, Lu^a and Co^b

Instructions for use can be found at www.blood.co.uk/reagents

	Rh	C	D	E	c	e	C ^w	M	N	S	s	P1	K	k	Kp ^a	Le ^a	Le ^b	Fy ^a	Fy ^b	Jk ^a	Jk ^b	Other	
1	R₁^wR₁	+	+	0	0	+	+	0	+	0	+	+	+	+	0	+	0	+	0	0	+	+	Cob+
2	R₂R₂	0	+	+	+	0	0	+	0	+	0	0	0	+	0	0	+	0	+	0	+	+	
3	rr	0	0	0	+	+	0	0	+	0	+	+	0	+	+	0	+	+	0	+	0	+	

EC REP Quality First International OÜ, Laki 30, 12915 Tallinn, Estonia