

### ID Panel Profile



IVD



Blood and Transplant

**NHSBT REAGENTS Panel 1**

Product	Lot No.	Product	Lot No.	Product	Lot No.
ID Panel in Alsevers	R144 3566	ID Panel Papainised in Alsevers	R154 3566	ID Panel in LISP	R146 3566
ID Panel in CellStab	R143 3566	ID Panel Papainised in CellStab	R153 3566	Expiry Date : 2024.08.08	
ID Panel in CellMedia	R163 3566	ID Panel Papainised in CellMedia	R173 3566		

Patient's Name	Ref. No.	Sample No.	Conclusion
D.O.B.			Tested by Date

Unless otherwise indicated, all cells are positive for Kp<sup>b</sup> and Lu<sup>b</sup> and negative for Wr<sup>a</sup> and Co<sup>b</sup>. \*Cell 9 is also Do(a-b+)

Instructions for use can be found at <http://www.blood.co.uk/reagents>

EC REP Quality First International OÜ, Laki 30, 12915 Tallinn, Estonia

	Rh	C	D	E	c	e	C <sup>w</sup>	M	N	S	s	P1	Lu <sup>a</sup>	K	k	Kp <sup>a</sup>	Le <sup>a</sup>	Le <sup>b</sup>	Fy <sup>a</sup>	Fy <sup>b</sup>	Jk <sup>a</sup>	Jk <sup>b</sup>	Other												
1	R <sub>1</sub> <sup>w</sup> R <sub>1</sub>	+	+	0	0	+	+	+	+	0	+	0	0	0	+	0	0	+	+	0	0	+													
2	R <sub>1</sub> R <sub>1</sub>	+	+	0	0	+	0	+	0	+	0	1	0	+	+	0	0	+	0	+	+	0	Do(a+b+)												
3	R <sub>2</sub> R <sub>2</sub>	0	+	+	+	0	0	0	+	0	+	1	0	0	+	0	+	0	+	0	+	0	HLA+ Bgb+												
4	r'r	+	0	0	+	+	0	+	0	+	0	2	0	0	+	0	0	+	+	0	+	+													
5	r''r	0	0	+	+	+	0	+	+	+	0	0	0	0	+	0	0	+	0	+	0	+													
6	rr	0	0	0	+	+	0	+	+	0	+	2	0	+	0	0	0	+	0	+	+	0													
7	rr	0	0	0	+	+	0	0	+	0	+	3	0	+	+	0	+	0	+	0	0	+	HLA+												
8	rr	0	0	0	+	+	0	+	0	0	+	2	0	0	+	+	+	0	0	+	+	0													
9	rr	0	0	0	+	+	0	0	+	+	0	1	0	0	+	0	0	+	0	+	0	+	* Cob+												
10	rr	0	0	0	+	+	0	+	0	+	0	0	+	0	+	0	0	0	0	+	0	+	Do(a-b+)												
Reagent																	DAT Profile	Anti-IgG	IgA	IgM	C3c	C3d	Ctrl												
Lot No.																	Result																		

(Template Version 01/04)