

### Papainised 3 Cell Screen Profile



**NHSBT REAGENTS**

|                               |                  |                   |
|-------------------------------|------------------|-------------------|
| Product                       | Lot No           | Expiry Date       |
| <b>Papainised in Alsevers</b> | <b>R124 3698</b> | <b>2024.08.29</b> |

Instructions for use can be found at [www.blood.co.uk/reagents](http://www.blood.co.uk/reagents)

|          | Rh  | C        | D        | E        | c        | e        | C <sup>w</sup> | P1       | K        | k        | Le <sup>a</sup> | Le <sup>b</sup> | Jk <sup>a</sup> | Jk <sup>b</sup> | Other |
|----------|---|----------|----------|----------|----------|----------|----------------|----------|----------|----------|-----------------|-----------------|-----------------|-----------------|-------|
| <b>1</b> | <b>R<sub>1</sub><sup>w</sup>R<sub>1</sub></b> | <b>+</b> | <b>+</b> | <b>0</b> | <b>0</b> | <b>+</b> | <b>+</b>       | <b>+</b> | <b>0</b> | <b>+</b> | <b>0</b>        | <b>+</b>        | <b>+</b>        | <b>+</b>        |       |
| <b>2</b> | <b>R<sub>2</sub>R<sub>2</sub></b>             | <b>0</b> | <b>+</b> | <b>+</b> | <b>+</b> | <b>0</b> | <b>0</b>       | <b>+</b> | <b>+</b> | <b>+</b> | <b>0</b>        | <b>+</b>        | <b>+</b>        | <b>0</b>        |       |
| <b>3</b> | <b>rr</b>                                     | <b>0</b> | <b>0</b> | <b>0</b> | <b>+</b> | <b>+</b> | <b>0</b>       | <b>+</b> | <b>0</b> | <b>+</b> | <b>+</b>        | <b>0</b>        | <b>0</b>        | <b>+</b>        |       |