

## rr Screen Profile Products PR106, PR107, PR108 & PR109



### NHSBT REAGENTS

Product	Lot No.	Expiry
<b>CellStab</b>	<b>R106 3283</b> <b>R107 3283</b>	<b>2024.11.14</b>
<b>CellMedia</b>	<b>R108 3283</b> <b>R109 3283</b>	

Unless otherwise indicated, all cells are positive for Kp<sup>b</sup> and Lu<sup>b</sup> and negative for Wr<sup>a</sup>, Lu<sup>a</sup> and Co<sup>b</sup>  
 Instructions for use can be found at [www.blood.co.uk/reagents](http://www.blood.co.uk/reagents)

		Rh	C	D	E	c	e	C <sup>w</sup>	M	N	S	s	P1	K	k	Kp <sup>a</sup>	Le <sup>a</sup>	Le <sup>b</sup>	Fy <sup>a</sup>	Fy <sup>b</sup>	Jk <sup>a</sup>	Jk <sup>b</sup>	Other	
<b>PR106/ PR108</b>	<b>1</b>	<b>rr</b>	<b>0</b>	<b>0</b>	<b>0</b>	<b>+</b>	<b>+</b>	<b>0</b>	<b>+</b>	<b>0</b>	<b>0</b>	<b>+</b>	<b>+</b>	<b>+</b>	<b>+</b>	<b>0</b>	<b>+</b>	<b>0</b>	<b>0</b>	<b>+</b>	<b>0</b>	<b>+</b>	<b>Cob+</b>	
	<b>2</b>	<b>rr</b>	<b>0</b>	<b>0</b>	<b>0</b>	<b>+</b>	<b>+</b>	<b>0</b>	<b>0</b>	<b>+</b>	<b>+</b>	<b>0</b>	<b>0</b>	<b>0</b>	<b>+</b>	<b>0</b>	<b>0</b>	<b>+</b>	<b>0</b>	<b>+</b>	<b>0</b>	<b>+</b>	<b>0</b>	
<b>PR107/ PR109</b>	<b>3</b>	<b>r'r</b>	<b>+</b>	<b>0</b>	<b>0</b>	<b>+</b>	<b>+</b>	<b>0</b>	<b>0</b>	<b>+</b>	<b>0</b>	<b>+</b>	<b>+</b>	<b>0</b>	<b>+</b>	<b>+</b>	<b>+</b>	<b>0</b>	<b>0</b>	<b>+</b>	<b>+</b>	<b>0</b>		
	<b>4</b>	<b>r''r</b>	<b>0</b>	<b>0</b>	<b>+</b>	<b>+</b>	<b>+</b>	<b>0</b>	<b>+</b>	<b>+</b>	<b>0</b>	<b>+</b>	<b>0</b>	<b>0</b>	<b>+</b>	<b>0</b>	<b>+</b>	<b>0</b>	<b>0</b>	<b>+</b>	<b>0</b>	<b>+</b>		

**EC REP** Quality First International OÜ, Laki 30, 12915 Tallinn, Estonia