

ID Panel Profile



IVD

NHSBT REAGENTS Panel 2



Blood and Transplant

Product	Lot No.	Product	Lot No.	Product	Lot No.
ID Panel in Alsevers	R141 3571	ID Panel Papainised in Alsevers	R15 n/a	ID Panel in LISp	R14 n/a
ID Panel in CellStab	R142 3571	ID Panel Papainised in CellStab	R152 3571	Expiry Date : 2024.12.19	
ID Panel in CellMedia	R162 3571	ID Panel Papainised in CellMedia	R172 3571		
Patient's Name D.O.B.		Ref. No.	Sample No.	Conclusion Tested by Date	

Unless otherwise indicated, all cells are positive for Kp^b and Lu^b and negative for Wr^a and Co^b.

Please note Cell 3 on the Papanised Cellmedia panel will be labelled as R172 5571.

Instructions for use can be found at <http://www.blood.co.uk/reagents>

EC REP Quality First International OÜ,Laki 30,12915 Tallinn, Estonia

	Rh	C	D	E	c	e	C ^w	M	N	S	s	P1	Lu ^a	K	k	Kp ^a	Le ^a	Le ^b	Fy ^a	Fy ^b	Jk ^a	Jk ^b	Other					
1	R ₁ ^w R ₁	+	+	0	0	+	+	0	+	+	0	0	0	0	+	0	0	+	+	0	+	0						
2	R ₁ R ₁	+	+	0	0	+	0	+	0	0	+	4	0	+	+	0	+	0	0	+	0	+						
3	R ₂ R ₂	0	+	+	+	0	0	+	0	+	0	3	0	0	+	+	0	+	0	+	+	0	HLA+					
4	r'r	+	0	0	+	+	0	0	+	0	+	0	0	0	+	0	+	0	+	+	+	+	0					
5	r''r	0	0	+	+	+	0	0	+	0	+	4	0	0	0	+	0	0	0	+	0	0	+					
6	rr	0	0	0	+	+	0	+	+	+	0	2	0	+	0	0	0	+	0	+	0	+						
7	rr	0	0	0	+	+	0	0	+	0	+	4	0	+	+	0	+	0	+	0	+	0	Cob+HLA+					
8	rr	0	0	0	+	+	0	+	0	+	0	0	0	0	+	0	0	+	+	0	+	0						
9	rr	0	0	0	+	+	0	0	+	+	0	0	0	0	+	0	+	0	0	+	0	+						
10	rr	0	0	0	+	+	0	0	+	0	+	1	+	0	+	0	0	+	0	+	+	0	Lub-					
Reagent																												
Lot No.																												
												DAT Profile		Anti-IgG		IgA		IgM		C3c		C3d		Ctrl				
												Result																

(Template Version 01/04)