

# ID Panel Profile



1434



Blood and Transplant

**NHSBT REAGENTS Panel 1**

Product	Lot No.	Product	Lot No.	Product	Lot No.
ID Panel in Alsevers	R144 3576	ID Panel Papainised in Alsevers	R154 3576	ID Panel in LISP	R146 3576
ID Panel in CellStab	R143 3576	ID Panel Papainised in CellStab	R153 3576	Expiry Date : 2024.12.26	
ID Panel in CellMedia	R163 3576	ID Panel Papainised in CellMedia	R173 3576		

Patient's Name		Ref. No.	Sample No.	Conclusion	
D.O.B.				Tested by	Date

Unless otherwise indicated, all cells are positive for Kp<sup>b</sup> and Lu<sup>b</sup> and negative for Wr<sup>a</sup> and Co<sup>b</sup>.

Instructions for use can be found at <http://www.blood.co.uk/reagents>

EC REP Quality First International OÜ, Laki 30, 12915 Tallinn, Estonia

	Rh	C	D	E	c	e	C <sup>w</sup>	M	N	S	s	P1	Lu <sup>a</sup>	K	k	Kp <sup>a</sup>	Le <sup>a</sup>	Le <sup>b</sup>	Fy <sup>a</sup>	Fy <sup>b</sup>	Jk <sup>a</sup>	Jk <sup>b</sup>	Other														
1	R <sub>1</sub> <sup>w</sup> R <sub>1</sub>	+	+	0	0	+	+	+	+	0	+	3	0	+	+	0	0	+	0	+	+	0															
2	R <sub>1</sub> R <sub>1</sub>	+	+	0	0	+	0	+	0	+	0	0	0	0	+	0	0	+	+	0	0	+															
3	R <sub>2</sub> R <sub>2</sub>	0	+	+	+	0	0	0	+	0	+	4	0	0	+	0	0	+	0	+	0	+															
4	r'r	+	0	0	+	+	0	+	0	0	+	0	0	0	+	0	0	+	0	+	+	0	HLA+														
5	r''r	0	0	+	+	+	0	+	+	+	0	2	0	0	+	0	0	+	+	0	+	0															
6	rr	0	0	0	+	+	0	0	+	0	+	4	0	+	0	0	+	0	0	+	+	+															
7	rr	0	0	0	+	+	0	+	0	+	0	2	0	+	+	0	+	0	+	0	0	+	Cob+														
8	rr	0	0	0	+	+	0	+	+	0	+	3	0	0	+	+	0	+	0	+	0	+															
9	rr	0	0	0	+	+	0	0	+	+	0	0	0	0	+	0	+	0	0	+	+	0	HLA+														
10	rr	0	0	0	+	+	0	0	+	0	+	4	+	0	+	0	+	0	+	0	0	+															
Reagent																	DAT Profile	Anti-IgG	IgA	IgM	C3c	C3d	Ctrl														
Lot No.																	Result																				