

### 3 Cell Screen Product Profile



#### NHSBT REAGENTS

Product	Lot No	Product	Lot No.	Product	Lot No.	Expiry Date
<b>Alsevers</b>	<b>R121 3705</b>	<b>CellStab</b>	<b>R122 3705</b>	<b>CellMedia</b>	<b>R123 3705</b>	2025.01.02

**Unless otherwise indicated, all cells are positive for Kp<sup>b</sup> and Lu<sup>b</sup> and negative for Wr<sup>a</sup>, Lu<sup>a</sup> and Co<sup>b</sup>**

Instructions for use can be found at [www.blood.co.uk/reagents](http://www.blood.co.uk/reagents)

	Rh	C	D	E	c	e	C <sup>w</sup>	M	N	S	s	P1	K	k	Kp <sup>a</sup>	Le <sup>a</sup>	Le <sup>b</sup>	Fy <sup>a</sup>	Fy <sup>b</sup>	Jk <sup>a</sup>	Jk <sup>b</sup>	Other
<b>1</b>	<b>R<sub>1</sub><sup>w</sup>R<sub>1</sub></b>	<b>+</b>	<b>+</b>	<b>0</b>	<b>0</b>	<b>+</b>	<b>+</b>	<b>+</b>	<b>0</b>	<b>+</b>	<b>+</b>	<b>+</b>	<b>+</b>	<b>+</b>	<b>0</b>	<b>0</b>	<b>+</b>	<b>0</b>	<b>+</b>	<b>+</b>	<b>0</b>	
<b>2</b>	<b>R<sub>2</sub>R<sub>2</sub></b>	<b>0</b>	<b>+</b>	<b>+</b>	<b>+</b>	<b>0</b>	<b>0</b>	<b>+</b>	<b>+</b>	<b>+</b>	<b>0</b>	<b>+</b>	<b>0</b>	<b>+</b>	<b>0</b>	<b>0</b>	<b>+</b>	<b>+</b>	<b>0</b>	<b>+</b>	<b>0</b>	
<b>3</b>	<b>rr</b>	<b>0</b>	<b>0</b>	<b>0</b>	<b>+</b>	<b>+</b>	<b>0</b>	<b>+</b>	<b>+</b>	<b>0</b>	<b>+</b>	<b>+</b>	<b>0</b>	<b>+</b>	<b>+</b>	<b>+</b>	<b>0</b>	<b>0</b>	<b>+</b>	<b>0</b>	<b>+</b>	

**EC REP** Quality First International OÜ, Laki 30, 12915 Tallinn, Estonia

(Template Version 01/04)