

### 3 Cell Screen Product Profile



#### NHSBT REAGENTS

Product	Lot No	Product	Lot No.	Product	Lot No.	Expiry Date
<b>Alsevers</b>	<b>R121 3707</b>	<b>CellStab</b>	<b>R122 3707</b>	<b>CellMedia</b>	<b>R123 3707</b>	2025.01.30

Unless otherwise indicated, all cells are positive for Kp<sup>b</sup> and Lu<sup>b</sup> and negative for Wr<sup>a</sup>, Lu<sup>a</sup> and Co<sup>b</sup>

Instructions for use can be found at [www.blood.co.uk/reagents](http://www.blood.co.uk/reagents)

	Rh	C	D	E	c	e	C <sup>w</sup>	M	N	S	s	P1	K	k	Kp <sup>a</sup>	Le <sup>a</sup>	Le <sup>b</sup>	Fy <sup>a</sup>	Fy <sup>b</sup>	Jk <sup>a</sup>	Jk <sup>b</sup>	Other
1	R <sub>1</sub> <sup>w</sup> R <sub>1</sub>	+	+	0	0	+	+	+	0	+	0	+	0	+	0	0	+	0	+	+	0	
2	R <sub>2</sub> R <sub>2</sub>	0	+	+	+	0	0	+	+	+	0	+	0	+	0	+	0	+	0	+	0	
3	rr	0	0	0	+	+	0	+	+	0	+	+	+	+	+	0	+	+	0	0	+	

EC REP Quality First International OÜ, Laki 30, 12915 Tallinn, Estonia