

3 Cell Screen Product Profile



NHSBT REAGENTS

Product	Lot No	Product	Lot No.	Product	Lot No.	Expiry Date
Alsevers	R121 3707	CellStab	R122 3707	CellMedia	R123 3707	2025.01.30

Please note that cell 2 in all diluents will be labelled as 4707

Unless otherwise indicated, all cells are positive for Kp^b and Lu^b and negative for Wr^a, Lu^a and Co^b

Instructions for use can be found at www.blood.co.uk/reagents

	Rh	C	D	E	c	e	C ^w	M	N	S	s	P1	K	k	Kp ^a	Le ^a	Le ^b	Fy ^a	Fy ^b	Jk ^a	Jk ^b	Other
1	R₁^wR₁	+	+	0	0	+	+	+	0	+	0	+	0	+	0	0	+	0	+	+	0	
2	R₂R₂	0	+	+	+	0	0	+	0	+	0	+	0	+	0	+	0	0	+	0	+	Cob+
3	rr	0	0	0	+	+	0	+	+	0	+	+	+	+	+	0	+	+	0	0	+	

EC REP Quality First International OÜ, Laki 30, 12915 Tallinn, Estonia

(Template Version 01/04)