

## 2 Cell Screen Profile PR101, PR102 & PR103



**NHSBT REAGENTS**

|                  |                  |                                |                  |
|------------------|------------------|--------------------------------|------------------|
| Product          | Lot No.          | Product                        | Lot No.          |
| <b>Alsevers</b>  | <b>R101 3552</b> | <b>CellStab</b>                | <b>R102 3552</b> |
| <b>CellMedia</b> | <b>R103 3552</b> | <b>EXPIRY DATE: 2025.04.03</b> |                  |

Unless otherwise indicated, all cells are positive for Kp<sup>b</sup> and Lu<sup>b</sup> and negative for Wr<sup>a</sup> , Lu<sup>a</sup> and Co<sup>b</sup>  
 Instructions for use can be found at [www.blood.co.uk/reagents](http://www.blood.co.uk/reagents)

|   | Rh                            | C | D | E | c | e | C <sup>w</sup> | M | N | S | s | P1 | K | k | Kp <sup>a</sup> | Le <sup>a</sup> | Le <sup>b</sup> | Fy <sup>a</sup> | Fy <sup>b</sup> | Jk <sup>a</sup> | Jk <sup>b</sup> | Other |
|---|-------------------------------|---|---|---|---|---|----------------|---|---|---|---|----|---|---|-----------------|-----------------|-----------------|-----------------|-----------------|-----------------|-----------------|-------|
| 1 | R <sub>1</sub> R <sub>1</sub> | + | + | 0 | 0 | + | 0              | + | 0 | 0 | + | +  | + | + | 0               | +               | 0               | 0               | +               | +               | 0               |       |
| 2 | R <sub>2</sub> R <sub>2</sub> | 0 | + | + | + | 0 | 0              | + | + | + | 0 | +  | 0 | + | 0               | 0               | +               | +               | 0               | 0               | +               |       |