

## rr Screen Profile Products PR106, PR107, PR108 &amp; PR109



Blood and Transplant

## NHSBT REAGENTS

| Product          | Lot No.                              | Expiry            |
|------------------|--------------------------------------|-------------------|
| <b>CellStab</b>  | <b>R106 3291</b><br><b>R107 3291</b> | <b>2025.06.26</b> |
| <b>CellMedia</b> | <b>R108 3291</b><br><b>R109 3291</b> |                   |

Unless otherwise indicated, all cells are positive for Kp<sup>b</sup> and Lu<sup>b</sup> and negative for Wr<sup>a</sup>, Lu<sup>a</sup> and Co<sup>b</sup>

Instructions for use can be found at [www.blood.co.uk/reagents](http://www.blood.co.uk/reagents)

|                         |          | Rh          | C        | D        | E        | c        | e        | C <sup>w</sup> | M        | N        | S        | s        | P1       | K        | k        | Kp <sup>a</sup> | Le <sup>a</sup> | Le <sup>b</sup> | Fy <sup>a</sup> | Fy <sup>b</sup> | Jk <sup>a</sup> | Jk <sup>b</sup> | Other       |
|-------------------------|----------|-------------|----------|----------|----------|----------|----------|----------------|----------|----------|----------|----------|----------|----------|----------|-----------------|-----------------|-----------------|-----------------|-----------------|-----------------|-----------------|-------------|
| <b>PR106/<br/>PR108</b> | <b>1</b> | <b>rr</b>   | <b>0</b> | <b>0</b> | <b>0</b> | <b>+</b> | <b>+</b> | <b>0</b>       | <b>+</b> | <b>0</b> | <b>+</b> | <b>0</b> | <b>+</b> | <b>+</b> | <b>+</b> | <b>0</b>        | <b>+</b>        | <b>0</b>        | <b>0</b>        | <b>+</b>        | <b>+</b>        | <b>0</b>        |             |
|                         | <b>2</b> | <b>rr</b>   | <b>0</b> | <b>0</b> | <b>0</b> | <b>+</b> | <b>+</b> | <b>0</b>       | <b>+</b> | <b>+</b> | <b>0</b> | <b>+</b> | <b>+</b> | <b>0</b> | <b>+</b> | <b>0</b>        | <b>0</b>        | <b>+</b>        | <b>+</b>        | <b>0</b>        | <b>0</b>        | <b>+</b>        |             |
| <b>PR107/<br/>PR109</b> | <b>3</b> | <b>r'r</b>  | <b>+</b> | <b>0</b> | <b>0</b> | <b>+</b> | <b>+</b> | <b>0</b>       | <b>0</b> | <b>+</b> | <b>0</b> | <b>+</b> | <b>+</b> | <b>0</b> | <b>+</b> | <b>0</b>        | <b>0</b>        | <b>+</b>        | <b>0</b>        | <b>w</b>        | <b>+</b>        | <b>0</b>        |             |
|                         | <b>4</b> | <b>r''r</b> | <b>0</b> | <b>0</b> | <b>+</b> | <b>+</b> | <b>+</b> | <b>0</b>       | <b>0</b> | <b>+</b> | <b>0</b> | <b>+</b> | <b>+</b> | <b>0</b> | <b>+</b> | <b>+</b>        | <b>0</b>        | <b>+</b>        | <b>+</b>        | <b>+</b>        | <b>0</b>        | <b>+</b>        | <b>Cob+</b> |

EC REP Quality First International OÜ, Laki 30, 12915 Tallinn, Estonia

(Template Version 01/04)