

Papainised 3 Cell Screen Profile



NHSBT REAGENTS

Product	Lot No	Expiry Date
Papainised in Alsevers	R124 3722	2025.07.31

Instructions for use can be found at www.blood.co.uk/reagents

	Rh	C	D	E	c	e	C ^w	P1	K	k	Le ^a	Le ^b	Jk ^a	Jk ^b	Other
1	R ₁ ^w R ₁	+	+	0	0	+	+	+	0	+	0	+	0	+	
2	R ₂ R ₂	0	+	+	+	0	0	+	0	+	+	0	+	0	
3	rr	0	0	0	+	+	0	+	+	+	+	0	w	+	