

Applicable Regulation	Product	LOT	2025.10.16
UK MDR (2002 as amended), IVDD	Reference Panel in Alsevers	R201 3319	
UK MDR (2002 as amended), IVDD	Reference Panel in CellStab	R202 3319	

Name	Hospital Number.	Conclusion	Tested by
D.O.B.	Sample No.		Date

	Rh	C	D	E	c	e	C ^w	M	N	S	s	P1	Lu ^a	Lu ^b	K	k	Kp ^a	Kp ^b	Le ^a	Le ^b	Fy ^a	Fy ^b	Jk ^a	Jk ^b	Wr ^a	Co ^b	Other					
1	R1wR1	+	+	0	0	+	+	+	0	+	0	0	0	+	0	+	0	+	0	+	+	0	+	0	0	0	0	Do(a+b-)				
2	R1wr'	+	+	0	0	+	+	0	+	0	+	0	0	+	0	+	0	+	+	0	0	+	0	+	0	0						
3	R1RZ	+	+	+	0	+	0	0	+	0	+	2	0	+	0	+	0	+	0	0	0	+	+	0	0	0	Do(a-b+)					
4	R1RZ	+	+	+	0	+	0	+	0	0	+	0	0	+	0	+	0	+	0	+	+	0	0	+	0	0						
5	Ror	0	+	0	+	+	0	+	0	+	+	3	0	+	0	+	0	+	0	0	0	0	+	0	0	0						
6	R1R1	+	+	0	0	+	0	0	+	0	+	4	0	+	+	0	0	+	0	+	+	0	0	+	0	0	HLA+					
7	R1R1	+	+	0	0	+	0	+	0	+	0	2	0	+	+	+	0	+	0	+	0	+	+	0	0	0						
8	R1R1	+	+	0	0	+	0	+	0	0	+	1	0	+	0	+	+	+	+	0	0	+	0	+	0	0						
9	R1R1	+	+	0	0	+	0	0	+	0	+	3	+	+	0	+	0	+	+	0	0	+	+	0	0	0						
10	R1R1	+	+	0	0	+	0	0	+	+	0	2	0	+	0	+	0	+	0	+	W	0	0	+	0	0						
																				DAT		Anti-IgG			IgA		IgM	C3c	C3d	Ctrl		



blood.co.uk/reagents



NHSBT Reagents, 14 Estuary Banks,
Speke, Liverpool, L24 8RB, Great Britain



Quality First International OÜ,
Laki 30, 12915 Tallinn, Estonia



UK
CA
0843



1434