

	<b>LOT</b>	
<b>CellStab</b>	R106 3304 R107 3304	2026.05.28
<b>CellMedia</b>	R108 3304 R109 3304	

		Rh	C	D	E	c	e	C <sup>w</sup>	M	N	S	s	P1	Lu <sup>a</sup>	Lu <sup>b</sup>	K	k	Kp <sup>a</sup>	Kp <sup>b</sup>	Le <sup>a</sup>	Le <sup>b</sup>	Fy <sup>a</sup>	Fy <sup>b</sup>	Jk <sup>a</sup>	Jk <sup>b</sup>	Wr <sup>a</sup>	Co <sup>b</sup>	Other
PR106 / PR108	1	rr	0	0	0	+	+	0	+	0	+	0	+	0	+	+	+	0	+	0	+	+	0	0	+	0	0	
	2	rr	0	0	0	+	+	0	+	+	0	+	0	0	+	0	+	0	+	+	0	0	+	+	0	0	0	
PR107 / PR109	3	r'r	+	0	0	+	+	0	+	+	+	+	+	0	+	0	+	0	+	0	+	+	+	+	0	0	0	
	4	r''r	0	0	+	+	+	0	+	+	0	+	+	0	+	+	+	+	+	0	+	+	+	+	0	0	0	



[blood.co.uk/reagents](http://blood.co.uk/reagents)



NHSBT Reagents, 14 Estuary Banks,  
Speke, Liverpool, L24 8RB, Great Britain



Quality First International OÜ,  
Laki 30,12915 Tallinn, Estonia



0843



1434

Controlled if copy number stated on document and issued by QA

(Template Version 03/02/2020)

Cross-Referenced in Primary Document: SOP883