

	<b>LOT</b>	
<b>CellStab</b>	R106 3305 R107 3305	2026.06.25
<b>CellMedia</b>	R108 3305 R109 3305	

		Rh	C	D	E	c	e	C <sup>w</sup>	M	N	S	s	P1	Lu <sup>a</sup>	Lu <sup>b</sup>	K	k	Kp <sup>a</sup>	Kp <sup>b</sup>	Le <sup>a</sup>	Le <sup>b</sup>	Fy <sup>a</sup>	Fy <sup>b</sup>	Jk <sup>a</sup>	Jk <sup>b</sup>	Wr <sup>a</sup>	Co <sup>b</sup>	Other
PR106 / PR108	1	rr	0	0	0	+	+	0	+	0	+	0	+	0	+	+	+	0	+	0	+	+	0	0	+	0	0	
	2	rr	0	0	0	+	+	0	+	+	0	+	0	0	+	0	+	0	+	+	0	0	+	+	0	0	0	
PR107 / PR109	3	r'r	+	0	0	+	+	0	+	0	+	+	+	0	+	0	+	+	+	0	+	0	+	+	0	0	0	
	4	r''r	0	0	+	+	+	0	+	+	+	0	+	0	+	0	+	0	+	+	0	0	+	0	+	0	0	



[blood.co.uk/reagents](http://blood.co.uk/reagents)



NHSBT Reagents, 14 Estuary Banks,  
Speke, Liverpool, L24 8RB, Great Britain



Quality First International OÜ,  
Laki 30,12915 Tallinn, Estonia



**UK  
CA**  
0843

**CE**  
1434

Controlled if copy number stated on document and issued by QA  
(Template Version 03/02/2020)

Cross-Referenced in Primary Document: SOP883