



NHS Blood & Transplant
500 North Bristol Park
Northway
Filton
Bristol
BS34 7QH

Antigen	CD233 / Band 3 (extracellular domain)
Clone	BRIC 200
Product Code	9468
Immunoglobulin Class	Mouse IgG1, kappa light chain

Tel: +44 (0)117 921 7200
Fax: +44 (0)117 912 5796
Web: <http://ibgri.blood.co.uk>

Antigen Description and Distribution

CD233 (also known as erythrocyte band 3, EPB3, anion exchange protein 1, AE1, solute carrier family 4, SLC4A1) is an integral membrane protein in human erythrocytes, present at approximately 10^6 copies per human erythrocyte. It comprises two domains that are structurally and functionally distinct. The cytoplasmic N-terminal 40kDa domain has binding sites for erythrocyte cytoskeletal proteins, namely ankyrin and protein 4.2, which help to maintain the mechanical properties and integrity of the erythrocyte. This domain also binds a number of other erythrocyte peripheral proteins. The 55kDa glycosylated C-terminal membrane-associated domain contains 12-14 membrane spanning segments which function as a chloride/bicarbonate anion exchanger involved in carbon dioxide transport. The cytoplasmic tail at the extreme C-terminus of the membrane domain binds carbonic anhydrase II. CD233 associates with the erythrocyte membrane protein glycophorin A (GPA) which promotes the correct folding and translocation during biosynthesis of CD233. Many CD233 mutations are known in man and these mutations can lead to two types of disease; destabilization of the erythrocyte membrane leading to hereditary spherocytosis, and defective kidney acid secretion leading to distal renal tubular acidosis. Other CD233 mutations that do not give rise to disease result in novel blood group antigens, which form the Diego blood group system. The CD233 gene is located on chromosome 17q21-q22¹.

Clone

BRIC 200 was made in response to intact human erythrocytes². BRIC 200 binds to an exofacial epitope on CD233 and agglutinates normal erythrocytes indirectly. BRIC 200 fails to agglutinate pronase-treated normal erythrocytes but agglutinates cells that have been treated with trypsin, sialidase, or 6% aminoethylisothiuronium bromide. The exact binding site of BRIC 200 is unknown but, because it fails to agglutinate pronase treated erythrocytes, it probably binds somewhere on the extracellular loop between amino acids 545 and 567². BRIC 200 gives markedly reduced titres with chymotrypsin treated cells. Chymotrypsin treatment of intact erythrocytes cleaves the membrane domain of CD233 after Try-553³ and Try-558⁴. BRIC 200 does not react with electrophoretically separated components of human erythrocyte membranes by immunoblotting. BRIC 200 has been used to study the coexpression of band 3 mutants and Rh polypeptides⁵. BRIC 200 has also been used to study the band 3 macrocomplex in the erythrocyte membrane⁶.

References

1. Bruce L.J. and Tanner M.J.A. (Review) *PROW* 2:9-17 (2001).
2. Smythe J.S., Spring F.A., Gardner B., Parsons S.F., Judson P.A. and Anstee D.J. *Blood* 85:2929-2936 (1995).
3. Steck T.L. *J Supramolec Struct* 8:311 (1978).
4. Jennings M.L., Adams-Lackey M., Denney G.H. *J Biol Chem* 259: 4652 (1984).
5. Beckmann R. *et al* (2001) *Blood* 97, 2496-2505.
6. Bruce L. J. *et al* (2003) *Blood* 101, 4180-4188.
7. Bruce L.J., Ring S.M., Anstee D.J., Reid M.E., Wilkinson S., Tanner M.J.A. *Blood* 85:541 (1995).
8. Groves J.D., Ring S.M., Schofield A.E., Tanner M.J.A. *FEBS Lett* 330: 196 (1983).
9. Toye A.M. *et al* (2004) *J. Cell Science* **117**, 1399-1410.